



Little **O**bjects **S**egmentation &
Tracking **E**nvironment

What is LOBSTER?

LOBSTER is an environment to batch process multidimensional microscopy images; it can identify biological objects and measure their morphology, dynamics and intensity distribution.

The software was designed with speed and simplicity in mind and a strong emphasis on data exploration and validation; it is fit for:

- Cell phenotyping in large high-content screening assays
- Object tracking in long microscopy time-lapses
- Filament tracing / spot counting in massive 3D images

Why LOBSTER?

A growing number of biological images analysis software are available, but most suffer from at least one of the following shortcomings:

- Too generic: no guidance to assemble algorithms into workflows
- Too specific: 2D images only, no object tracking, etc.
- Not scalable: limited image size, slow, too sensitive to varying quality
- Not flexible: fixed workflow, difficult to customize to a specific project
- Not expressive: too much code to write, prone to errors


Highlights

- ✓ Handle multi-dimensional images (3D + channels + time)
- ✓ Virtually no image size limit
- ✓ Consistent set of robust, documented IA algorithms (functions)
- ✓ Easily assemble functions into readable IA workflows (journals)
- ✓ Many sample applications (images + journals)
- ✓ Image analysis server: queue jobs, email results

Usability


- ✓ Simple to deploy (Matlab), no programming (plumbing only)
- ✓ Streamlined flow, only 4 basic type of image masks
- ✓ Open, layered architecture, human readable data at each step
- ✓ Detailed documentation built as tutorials
- ✓ Images, results and workflows accessible from hyperlinks

Architecture








Sample Applications



Object segmentation



Object tracking




Filament / tube tracing




3D Object detection




Spot detection




Filament segmentation





**Machine learning
texture classification**





**Machine learning
nuclei detection**





Cell detection





Out of focus resilient object segmentation








Clustered 3D object segmentation




Surface segmentation



Cell detection over 3D mesh







Objects extraction from seeds





Robust convex particles splitting

Advanced denoising





Membrane restoration





Destriping

Restoration





**Phase contrast
restoration**


Deconvolution





**Tube to filament
restoration**





Extended depth of view






Speed
measurement
(PIV)



Orientation measurement



Points clustering

Image Features Measurements

Calling Journals: JENI

(Journal ENgine Interface)

```
>> [InputFolder MaskFolder] = JENI ('NucleiCytoo_GradWaterTilesMerge.jl');
```

Input and output image folders (MATLAB variables)

Journal to call

```
Journal: E:\LOBSTER/Journals/jl/NucleiCytoo_GradWaterTilesMerge.jl-->Launch
Input Folder: ./Images/NucleiCytoo/
Output Folder: ./Results/Images/NucleiCytoo/
Warning: Output folder not empty!!
> In JENI_Images (line 53)
  In JENI (line 95)

-----
Found 2 image(s) / group(s)
Rescale: 1
MinLocalFocus: 0
ExportDist: 0

-----
processing... 100% [.....]
Time spent inside function fxg_mGradWaterTiles : 0.11 s
Time spent inside function fxm_lTilesMerge : 0.39 s
Total processing time: 0.51 s
```

Open journal

Rerun journal

Parsed files
+ configuration







Execution time



Open folder



Open folder

Open file

Open file

Name	Date
 NucleiCytoo1_C00.tif	06/05/2015 13:24
 NucleiCytoo1_C01.tif	06/08/2017 21:46
 NucleiCytoo1_C02.tif	04/12/2017 18:30
 NucleiCytoo3_C00.tif	31/01/2017 18:22
 NucleiCytoo3_C01.tif	06/08/2017 21:46
 NucleiCytoo3_C02.tif	04/12/2017 18:31

Name	Date
 NucleiCytoo1_C00_L2.tif	26/02/2018 14:24
 NucleiCytoo3_C00_L2.tif	26/02/2018 14:24



Measurements: IRMA

(Image Regions Measurements & Analysis)

Coming from JENI

```
>> ReportFolder = IRMA (MaskFolder, '.', 'Objs', 2, 1, InputFolder, '*C01*.tif');
```

Report folder

- Mask images to measure
- File filter

- Measurement type
- 2D
- Zratio = 1

Intensity channels (optional)

```
Intensity measurements in 1 channels
Intensity channels folder: ./Images/NucleiCytoo/
Found 2 datasets in input folder
```

```
-----
Processing dataset 1 of 2
Number of objs:      58      Sum area      : 24568
Mean intensity:     103.71    Non null pix: 24568
```

```
-----
Processing dataset 2 of 2
Number of objs:       2      Sum area      : 781
Mean intensity:      75.00    Non null pix: 781
```



```
-----
Masks folder: ./Results/Images/NucleiCytoo/
Chan folder:  ./Images/NucleiCytoo/
Reports folder: ./Results/Reports/NucleiCytoo/
```

Parsed files

Measurements
summary

	Area	Centroid_1	Centroid_2	BoundingBox_1	BoundingBox_2	BoundingBox_3	BoundingBox_4	Meanint	NonNullPix
1	375	30.50	26.87	20.50	14.50	22	26	50	375
2	485	33.56	126.00	20.50	112.50	25	26	50	485
3	410	32.57	301.09	21.50	289.50	23	23	50	410
4	423	36.46	69.86	22.50	58.50	30	23	50	423
5	485	36.56	204.00	23.50	190.50	25	26	50	485
6	492	43.76	154.51	30.50	142.50	25	25	50	492
7	492	46.76	232.51	33.50	220.50	25	25	50	492
8	512	52.65	82.52	35.50	68.50	32	26	50	512

Open folder

Name	Date modified
 NucleiCytoo1_C00_L2.csv	01/03/2018 16:37
 NucleiCytoo3_C00_L2.csv	01/03/2018 16:37

Open file

Exportation: JOSE



(JOB Scene Exporter)

Coming from JENI Coming from IRMA

```
>>JOSE(InputFolder,'*_C00*',ReportFolder,'Objs',ReportFolder,'IJ',);
```

Link C00 from original images Link Objects bounding boxes Export scene to this folder for IJ (SceneViewer macro)

.sce files parsed by IJ
SceneViewer macro
from report folder




3 steps together: job file (*.m)

```
[InputFolder MaskFolder] = GENI('NucleiCytoo_GradWaterTilesMerge.j1');  
ReportFolder = IRMA(MaskFolder, '.', 'Objs', 2);  
JOSE(InputFolder, '*_C00*', ReportFolder, 'Objs', ReportFolder, 'IJ');
```


Job 'NucleiCytoo_scene.m'

Journals


```
InputFolder = './Images/NucleiCytoo/';  
OutputFolder = './Results/Images/NucleiCytoo/';  
Fill = 1; display filled masks  
Lbl = 1; display labeled masks (colors)  
  
@iA = '*C00*.tif'; image filter (InputFolder)  
  
@fxg_mGradWaterTiles [iA] > [L];  
params.GaussianRadInt = 2;  
params.ExtendedMinThr = 2;  
/endf first function  
  
@fxm_lTilesMerge [L, iA] > [L2];  
params.GaussianRad = 2;  
params.MinObjArea = 175;  
params.MinSal = -0.5;  
params.MaxValleyiness = 1.075;  
params.ConcavityThresh = 0.4;  
/endf second function  
  
/show iA > L2; display iA + L2 overlay (image viewer)  
/keep L2 > tif; export L2 to tif (OutputFolder)
```



iA



L



L2

Journal 'NucleiCytoo_GradWaterTilesMerge.jl'

Documented Functions

```
>> help fxm_lTileMerge
```

```
Merge/split particles in binary mask to reconstruct loosely convex bright objects.  
The algorithm also requires the original image to use intensity information.
```

```
Sample journal: NucleiCytoo\_GradWaterTilesMerge.jl
```

← Hyperlink to sample journal
demoing the use of this function


```
Input: 2D binary mask, 2D original image  
Output: 2D label mask
```

```
Parameters:
```


GaussianRad:	Gaussian blur pre-filter radius (pix)
MinObjArea:	Minimum object area (pix)
MinSal:	Minimum tiles saliency (low -0.5 -> 0.5 high)
MaxValleyiness:	Maximum tiles valleyiness (rescue more 0.75 -> 1.25 rescue less)
ConcavityThresh:	Concavity threshold (sensitive 0.25 -> 0.5 coarse)

Basic Measurements


Object mask (binary)



Object label mask



Seed mask



Skeleton mask



Objects

Area
Centroids (X,Y,Z)
Bounding boxes (6 coords)
Mean intensity

Objects

Area
Centroids (X,Y,Z)
Mean intensity




Seeds

Centroids (X,Y,Z)
Mean intensity

Skeletons




Mask volume
Length
Nb brch/end points
Object volume
Image volume
Mean intensity

Advanced Measurements



Filaments:
Branch length distribution
Width distribution

Secondary objects count **Objects distance distribution**




Objects colocalization **Statistical colocalization** **Clustering**


- Normalized objects overlap
- Seeds FP / FN / TP rates
- Predicted ++ probability (assuming random seeds)
- Actual ++ observation prob. (assuming random seeds)
- Colocalization evidence
 $P(r+ \mid g+) / P(r+ \mid g-)$
- Exclusion evidence
 $P(r- \mid g+) / P(r- \mid g-)$
- Clustering factors (r and g)

Exportation & validation

... in massive images




**Display object bounding
Box + color-code (IJ)**




**Montager IJ macro:
Auto-montage objects
+ link to their locations**


3D models exportation




Seeds (xls)
CellInsight (IJ)



Filaments (SWC)
Neuromantics



Tracks (SWC)
Neuromantics



Objects (STL)
MeshLab


Image Analysis server: JULI

(Job UpLink)


JULI enabled job

```
>> JULI (PathToMonitoredFolder, PathToLogFile, 'user@gmail.com', 'password');
```


```
InputFolder = 'E:/LOBSTER_sandbox/Images/FISH/';  
Dstmail = 'youremail@youremailprovider.com';  
[InputFolder MaskFolder1] = GENI('FISH_nucseg.jl', InputFolder, 1);  
[InputFolder MaskFolder2] = GENI('FISH_sptdet.jl', InputFolder, 2);  
AttachmentFolder = IRMA(MaskFolder1, 1, 'Objs', 2, 1, MaskFolder2);
```



Monitored
Folder (FTP)




LOBSTER + FTP + mail server
(MATLAB + Filezilla server)



Attachment Folder (csv, scenes, models)
zipped and sent to Dstmail

OUTLOOK



- 3D annotations edition tool (partly done)
- Deep learning all rounder module: images -> masks (in works)
- Job builder: graphical wizard
- Whatever else you find useful...

--> sebastien.tosi@irbbarcelona.org